

The Large Plasma Chamber

John Meyer

Deutsches Zentrum für Luft- und Raumfahrt
Weßling, Germany

Special thanks to:
Vladimir Nosenko, Hubertus Thomas




Knowledge for Tomorrow




Large Chamber Setup at DLR




Setup Diagram





Experimental Parameters



Polymer (MF) microspheres:


Diameter	$d \approx 9 \mu\text{m}$
Charge	$Q \approx -2 \cdot 10^4 e$
Spacing	$a \approx 0.9\text{-}1.6 \text{ mm}$
Number	$\approx 14,000$

Argon RF plasma:


pressure	2.50-5.00 Pa
RF power	100-420 W
Self Bias	-10- -40 V




Electrode Current and Voltage






Chamber Power vs Input




Self Bias











2.75 Pa
200 W



Particle Spacing vs Power**Particle Spacing vs Pressure**







5 Pa 125 W



5 Pa 175 W



Sample Spectra at Powers



Longitudinal

Transverse

Together




Longitudinal

Transverse

Together






Sample Spectra at Pressures



Longitudinal

Transverse

Together




Longitudinal

Transverse


Together





Q vs Power



Q vs Pressure



Screening Length vs Power**Screening Length vs Pressure**

Future Work

- Langmuir Probe measurements
- Bulk behavior analysis
- System-size dependence of thermal conductivity
- System-size dependence of dispersion relations
- Mode-coupling instability
- Possible photophoretic force experiment



Bulk Behavior




FIG. 1. Macroscopic shielding parameter $\bar{\kappa}=a/\lambda$ vs the Debye shielding parameter $\kappa=r_0/\lambda$ [Eq. (3)] for the finite Yukawa disk with $n=100, 180, 320, 560, 1000, 1800$, and 3200 particles.




FIG. 2. Percentage error between the squared breathing frequency computed using Eq. (5) and exact eigenvalues vs the number of particles n and the Debye shielding parameter κ .




FIG. 3. Normalized, squared breathing frequency vs the number of particles n computed using the exact particle model for values of the Debye shielding parameter $\kappa=0.2$ to 10 .

T. E. Sheridan, POP 14, 032108 (2007)



Thank You

