

Remote Raman detection of chlorine with deep UV excitation wavelengths

Emanuela C. A. Gallo, Frank Duschek

German Aerospace Center, Institute of Technical Physics



2nd Scientific International Conference on CBRNe

SICC Series Conference

Rome – Italy, December 10-12 2020



Knowledge for Tomorrow



Motivation

- Remote measurements to detect chlorine gas
- Chlorine gas exposure:




- **Industrial**
manufacturing
storage
handling



Motivation

- Chlorine gas exposure:

➤ **Warfare agent** - chlorine, mustard gas, bromine and phosgene



- Lethal doses:

Chlorine	6000 [mg min /m ³] - 2070 ppm min
Mustard gas	900 [mg min /m ³] - 230 ppm min
Tear gasses	30 mins temporary effect
Phosgene	3000 [mg min /m ³] - 740 ppm min




Introduction

- Remote detection set up optimized (distance of 60 cm)
- Change excitation WL in the deep UV => To maximize Cl₂ signal
- Detection limits in acquisition times
- Solution for background interference => Lab test: Cl₂ must be enclosed
- **SYSTEM SETUP**
- **SAMPLE SETUP**
- **RESULTS**
- **CONCLUSION AND FUTURE DEVELOPMENTS**



System setup




- Nd:YAG laser
- dye laser
- mixing unit
- 190-900 nm, 2.5 mJ/pulse, 10Hz
- Liq. N₂ cooled spectr. 2400 grooves/mm
- 60 cm remote distance detection
- Laser filter
- Cl₂ ~ 1 mg/ml




Sample setup

- Self-made system to fill a chlorine gas cell
- chlorine pure gas tank 99.8%
- Vacuum pump below 0.4 bar
- Return line (chemical dechlorinators) sodium thiosulfate, sodium hydroxide, deionized water
- local $\Delta P < 0.5$ bar to avoid window breaking limits
- final sample concentration was 0.36 bar




Results

Absorbance vs wavelength



- Edinburgh Instruments FS5 Spectrofluorometer
- 1nm step, 0.5 s sample time per wavelength
- max at 330 nm, FWHM 60 nm
- molar extinction coefficient, path length:
[Cl₂] = 0.36 bar

Fluorescence vs wavelength




- standard UV laser sources
- negligible fluorescence



Results

Background signal




empty quartz glass cell
as ref. background

- laser energy density below 20 mJ/cm²
- 224, 232, 235 nm tested to maximize signal
- broad peak 410 cm⁻¹ + sharp one 490 cm⁻¹ , in agreement with literature




Results

background & Cl₂ signal



Cl₂ signal




- quartz cell material overlapping with chlorine Raman signal
- Cl₂ sharp peak at 554 cm⁻¹ (15 cm⁻¹ FWHM)
- expected for Cl₂ at 554, 547, and 539 cm⁻¹
- broad peak quartz residual, in agreement with literature
- hard to separate the two
- remove cell not possible => change material




Results

- Avoid cell material interference around region of interest => Raman grade CaF_2 , diamond



- Signal increases at higher excitation wavelengths for both cell material and Cl_2
- strong overlapping of the unwanted quartz material with chlorine at higher wavelengths



- Cl_2 growing increasing the incoming laser wavelength
- Increasing time increases signal by 24%
- 50 s detection limit (peak interference, not intensified camera)



Conclusions & Future Developments

- Chlorine gas was detected in a remote Raman configuration: not standard setup
- Collimated configuration, 60 cm detection distance,
Total cell 1 mg/ml, probed volume (125 μg)
laser energy density < 20 mJ/cm^2 , acquisition time 50 s
- Reduce acquisition time using an ICCD detector
- Avoid strong quartz interference around region of interest (peaks overlapping) => Raman grade CaF_2 , diamond
New windows material required
- excitation wavelengths, concentration limits



Thank you!

Questions



Remote Raman detection of chlorine with deep UV excitation wavelengths

Emanuela C. A. Gallo, Frank Duschek

German Aerospace Center, Institute of Technical Physics

email: emanuela.gallo@dlr.de



2nd Scientific International Conference on CBRNe

SICC Series Conference

Rome – Italy, December 10-12 2020



Knowledge for Tomorrow