

# UV Raman Spectroscopy for Explosives Detection

E. Gallo, C. Kölbl, A. Walter, A. Köhntopp, F. Duschek  
W. Schweikert, F. Schnürer



DLR Lampoldshausen, Institute of Technical Physics  
Fraunhofer Institute for Chemical Technology

## ➤ Objective:

- Identification of explosive traces at a safe distance using UV Raman spectroscopy
  - Determination of the detection limits for samples on glass substrates.  $\mu\text{g}$  amounts resemble a chemical trace left by a fingerprint


## ➤ Lab setup:

- Portable 355 nm, 224 nm lasers
  - Collecting optics: UV lens ~ 60 cm in front of the sample - spectrometer (liquid  $\text{N}_2$  cooled CCD, grating 2400 grooves/mm, laser line filter)
  - Photo degradation minimized (laser energy <  $10 \text{ mW}/\text{cm}^2$ )
- UV Laser beam features:  
1 mW, 3 s minimum acquisition time




## ➤ Sample preparation:

- Standardized samples of milligrams and traces ( $\mu\text{g}$  range) were ink-jet printed on a glass substrate
- Samples simulate possible explosive trace amount adhering on a surface





## ➤ Results:





— RDX acetonitrile sol 1mg/ml  
— PETN acetonitrile sol 1mg/ml  
— Laser WL 224nm  
— Laser WL 355nm

- PETN (red) and RDX (blue), samples <  $100 \mu\text{g}$ , 20 s, 355 nm:



- Minimized acquisition time 3 s, 355 nm, <  $100 \mu\text{g}$  and <  $50 \mu\text{g}$ :



## ➤ Findings:

- A UV Raman setup for measuring explosive traces was successfully tested
- Lowest detection limit was <  $50 \mu\text{g}$  for 0.6 m distance for 3 s, at 355 nm
- No distinguishable signal for laser at 224 nm: too close to absorption peaks
- Tested compounds successfully identified

## ➤ Future developments:

- Find the optimum excitation wavelength for later applications (i.e. miniaturization)
- Test different surfaces, inhomogeneous contamination, low concentrations
- Eye safe range



Deutsches Zentrum  
für Luft- und Raumfahrt  
German Aerospace Center

Corresponding author:  
[emanuela.gallo@dlr.de](mailto:emanuela.gallo@dlr.de) +49 6298 28 494



Fraunhofer



Knowledge for tomorrow