

3 μm water vapor self- and foreign-continuum: New method for determination and new insights into the self-continuum

M. Birk, G. Wagner, J. Loos, K.P. Shine



Knowledge for Tomorrow



Determination of H₂O continua in the past

Method

- Transmittance spectra from FTS, absorption coefficients from CRDS
- Calculation of H₂O local lines using mostly HITRAN parameters
- Only use troughs, subtract local line contributions

Disadvantages (especially for FC)

- Large gaps
- Low resolution
- Errors due to inadequacy of line parameter database



New method for H₂O continuum determination

Method

- Usage of same spectra for line parameter and continuum determination
- Continuum extracted from fitted baseline from microwindow-based multispectrum fit

Advantage

- Continuum information not constrained to troughs – only few gaps
- Local line errors small
- Sufficient spectral resolution – determined by microwindow width, baseline = polynomial, degree ≤ 2
- Advantage is largest for in-band FC

Prototyping

- Water measurements between 1800 and 4000 cm⁻¹ with high quality line parameters already determined [Loos2017;Loos2017a]





- Since 1990 high resolution spectroscopy for spectroscopic databases of atmospheric constituents
- Commercial Bruker IFS125 HR FT spectrometer 10-40000 cm^{-1}



Multireflection cell

- 14-165 m absorption path, accuracy 0.1%
- 190-350 K temperature range
- Mirrors actively thermalized
- 0.1 K temperature homogeneity
- Attached to Bruker IFS 125 HR
- Transfer optics in sample chamber
- Mirror coating selected for minimum reflection loss



CH₄
H₂O
CO CO₂





Other infrastructure: Gas handling

- Precision pressure transducers: mks Baratron 0.05% accuracy
- Temperature sensors: Lakeshore Pt100, 0.1 K accuracy
- 800 l stainless steel gas mixing chamber, stirrer
- Flow and pressure controllers



Multispectrum fitting software



Line Explorer window showing preferences and fit results:

GLOBAL PREFERENCES #1

Select spectra: ☒ 0 ☒ 1 ☒ 2 ☒ 3 ☒ 4 ☒ 5 ☒ 6 ☒ 7 ☒ 8 ☒ 9 ☒ 10 ☒ 11

SPECTRA PREFERENCES

of baseline parameters: 0

☐ Fit Channelling A: 0.0000000 T: 1.0000000 maxWN: 1606.7500

☐ fit cal factor

ILS: sincBox

GLOBAL PREFERENCES #2

Fit rises (only when sincBox is selected)

☐ 1.2650000

1.2650000

Start Fit

fit OK

Residue Factor: x10

Plot Absorption

LINE PREFERENCES

0

11 1606.716310 1.036E-21 9.408E-02 0.9180.450 447.25240.720.001040

line profile: Speed-dependent Voigt

select parameters for SDV fit of line: 0

<input checked="" type="checkbox"/> sigma	1606.7163	5.0552128e-006	0.0000 %
<input type="checkbox"/> S	1.0360000e-021	0.00000000	0.0000 %
<input checked="" type="checkbox"/> gamma0_air	0.091812379	4.5797314e-005	0.0499 %
<input checked="" type="checkbox"/> gamma2_air	0.013911432	0.00018493259	1.3291 %
<input type="checkbox"/> eta0_air	0.001040000	0.00000000	0.0000 %
<input type="checkbox"/> eta2_air	0.00000000	0.00000000	NaN %
<input type="checkbox"/> delta0_air	0.00000000	0.00000000	NaN %
<input type="checkbox"/> n_air	0.72000000	0.00000000	0.0000 %
<input type="checkbox"/> gamma0_self	0.45000000	0.00000000	0.0000 %
<input type="checkbox"/> gamma2_self	0.00000000	0.00000000	NaN %
<input type="checkbox"/> eta0_self	0.001040000	0.00000000	0.0000 %
<input type="checkbox"/> eta2_self	0.00000000	0.00000000	NaN %
<input type="checkbox"/> delta0_self	0.00000000	0.00000000	NaN %
<input type="checkbox"/> n_self	0.72000000	0.00000000	0.0000 %

Automatic parameter selection SET Initial Pars

☐ equal gamma2


☐ equal eta0

☐ equal eta2


☐ equal delta0

☐ equal n

Multispectrum fit of pure H2O spectra



Multispectrum fit of air-broadened H₂O spectra



Measurements used for continuum determination


- $SC \propto p_{H_2O}^2 \cdot l = IC_s$
- $FC \propto p_{H_2O} \cdot p_{air} \cdot l = IC_f$

#	p_{H_2O}/mbar	p_{tot}/mbar	l/m	IC_s	IC_f
1	4.9626	4.9626	14.50	350	
2	20.0	20.0	14.50	5800	
3	4.9651	4.9651	72.08	1773	
4	20.0	20.0	72.08	28800	
5	4.9655	4.9655	177.3	4354	
6	20.0	20.0	176.7	70680	
7	6.0	300	72.08	2592	129600
8	2.0	1000	72.08	288	144000
9	7.5	1000	72.08	4050	540000
10	6.0	100	168.21	6048	100800
11	6.0	300	168.21	6048	302400
12	4.8	800	168.21	3870	645000




Transmittance contributions of pure H₂O spectra

green: observed, red: pure line spectrum, black: baseline




Transmittance contributions of air-broadened H₂O spectra

green: observed, red: pure line spectrum, black: baseline




Self-Continuum multispectrum fit to baselines

black: observed, red: calc, green: OMCx10




Foreign-Continuum multispectrum fit to baselines

black: observed, red: calc, blue: SC contribution, yellow: FC contribution, green: OMCx10




Validation of SC with air-broadened measurements




Comparison of SC with CAVIAR

lower traces: continua, upper traces: uncertainties/differences x 4,
zero blue line. Black: present work, green: CAVIAR, red: CAVIAR –
present work



Comparison of FC with CAVIAR and MT_CKD3.2




Shape of self-continuum

- H₂O-dimer dissociating into monomer
- Spectrum in between bound dimer and monomer
- Ab initio can calculate spectra but high energy level density close to dissociation limit
- Simpler approach: Bound dimer + quasibound dimer
- Quasibound dimer has shape of monomer band but large broadening due to short lifetime (FWHM=14-40 cm⁻¹, Ptashnik, 2011).
- New concept for estimating bound dimer spectrum: Generic band shapes of parallel and perpendicular rovibrational bands



Shape of bound dimer bands



- Strongest vibrational fundamentals in 3 μm region from ab initio/experiments:
 - symmetric OH stretch of the donor H_2O : ν_3
 - antisymmetric OH stretch of the donor H_2O : ν_1
 - antisymmetric OH stretch of the acceptor H_2O : ν_9
- Ground state rotational constants: $A \approx 7 \text{ cm}^{-1}$, $B \approx C \approx 0.2 \text{ cm}^{-1}$



Mukhopadhyay, 2015


Shape of bound dimer bands

- Ptashnik: All bands have Lorentzian shape with 60 cm^{-1} FWHM
- New approach: Generic band shapes of parallel and perpendicular bands of symmetric top
- ν_3 , ν_1 : dipole moment changes parallel and perpendicular to A-axis, ν_9 : perp.
- Effects smearing band structure: diff. rot. const. in ground and excited vib. states, centrifugal distortion, hot bands (ca. 85% of total intensity), asymmetry splitting, torsional tunneling splitting, lifetime broadening, Cave: floppy molecule
- FWHM parallel band rovibrational lines: 5 cm^{-1} , perpendicular band 10 cm^{-1}



Fit of SC continuum

#	K_{Db} /atm ⁻¹	K_{Dq} /atm ⁻¹	K_{eq} /atm ⁻¹	$f_{S_{v3}}$	$f_{S_{v1}}$	$f_{S_{v9}}$	$f_{para_{v3}}$	$f_{para_{v1}}$	$\nu_{0,v3}$ /cm ⁻¹	$\nu_{0,v1}$ /cm ⁻¹	$\nu_{0,v9}$ /cm ⁻¹	FWHM qD /cm ⁻¹
fit	0.025	0.042	0.067	0.48	0.27	0.25	0.54	0.6	3610	3724	3749	20
lit			0.052	0.55(8)	0.28(6)	0.18(3)			3597	3730	3749	



Conclusion

- New method for H₂O in-band continuum determination prototyped, using same spectra for line parameters and continuum
- Example: FC and SC in 3 μm region
- Continua show better spectral coverage and resolution, and smoothness
- SC in good agreement with CAVIAR work
- FC in agreement with CAVIAR work, showing P-, Q-, and R-branch
- SC was fitted applying bound and quasibound dimer contributions
- Rovibrational bands of bound dimer were simulated by generic parallel and perpendicular symmetric top bands
- Two double peaks in SC can be modelled with parallel rovibrational bands
- Dimer band positions and strengths are in reasonable agreement with theoretical work

