

Multispectrum fit of non-Voigt lineshape in the H₂O v₂ band

J. LOOS^a, M. BIRK^a, G. WAGNER^a, F. HASE^b, D. DUBRAVICA^b, M. PALM^c, A. SADEGHI^c

^a Deutsches Zentrum für Luft- und Raumfahrt (DLR), Institut für Methodik der Fernerkundung, joep.loos@dlr.de

^b Karlsruher Institut für Technologie (KIT), Institut für Meteorologie und Klimaforschung

^c Universität Bremen, Institut für Umweltphysik (IUP)

Abstract

A new fitting tool for analysis of multiple molecular absorption spectra utilizing a microwindow-based line-by-line-approach has been developed. Its capabilities include the choice of numerous different line shape models, from a simple Voigt to more sophisticated models like a speed-dependent Galatry including line-mixing. A comfortable manual mode as well as a fully automatic mode have been implemented including various quality assessment procedures like the monitoring of correlation coefficients or the supply of useful information e.g. needed for file cuts (single spectrum residuals) [1].

As a first application the new tool is used to re-analyze water vapor absorption spectra in the 1250-1750 cm⁻¹ range [1,2]. The measurements include pure water as well as water/air-mixture measurements and cover a wide range of column densities. The total air pressure and partial pressure ranges were 50-1000 mb and 0.001-5 mb, respectively.

Whereas the original analysis was based on single spectrum fits applying the Voigt procedure, in the present multispectrum fit the speed-dependent Voigt lineshape was used. The advantages of a multispectral analysis approach as well as the need for consideration of narrowing effects is illustrated by the presentation of differences in residuals as well as resulting line parameters for selected transitions. As indicated in [3] opaque as well as non opaque lines could be fitted with the speed-dependent Voigt while the pure Voigt yields to narrow opaque lines.

New multispectrum fitting tool

The software - written in IDL - has been developed for fitting of multiple absorption spectra recorded with a Fourier-transform spectrometer. Great effort was taken to make it as generic and comfortable as possible. Although it is tailored to the needs of high resolution FT-spectroscopy, it might be easily customized to be applied to spectra recorded with other instrument types. The results of the tool were validated vs. a single-spectrum IDL DLR fitting tool [4] and FITMAS [5].

Line models

- Voigt
- Speed-dependent Voigt [6]
- Speed-dependent Galatry + Rosenkranz line mixing (based on PROFFWD-routine [7,8])

Spectrum/microwindow-specific fitting capabilities

- polynomial baseline
- channelling
- wavenumber-calibration factor
- offset
- ILS (delta, sincbox, modulation + phase as a function of path difference as provided by LINEFIT [9])

Multispectrum fitting capabilities

- line parameters (position, intensity, self-/foreign-broadening + T-dependency, self-/foreign-shift + T-dependency, self-/foreign-narrowing parameters, line-mixing Y-parameters)
- entrance aperture

Automatic mode

- automatic microwindow-, spectra- and line-fit parameter selection
- automatic fit
- automatic iteration of fit parameters (e.g. parameter error above threshold)

Quality assessment

- goodness-of-fit (reduced χ^2 , residua analysis planned)
- identification/prevention of parameter correlations

$$C_{ij} = \frac{(J^T W J)^{-1}_{ij}}{\sqrt{(J^T W J)^{-1}_{ii} (J^T W J)^{-1}_{jj}}} \quad |C_{ij}| \leq C_{\max}$$

i, j = fitted parameter number

- single-spectrum-fits with line-specific parameters (position, intensity, width, narrowing)
- differences of calculated parameters from multispectrum fit results
- measurement-specific differences can be investigated (as a function of position, intensity, ...) -> file cut [1], e.g. Fig. 1
- helps to identify systematic spectrum-specific errors
- calculation of contribution matrix B
- percentage information contribution of single spectrum to parameter

$$B_{sp}^i = \frac{\sum_k [(J^T W J)^{-1} J^T W]_{kp}}{\max_k [(J^T W J)^{-1} J^T W]_{kp}} \quad B_{sp} = \sum_s B_{sp}^i$$

W = weight matrix, s = spectrum number,
 p = parameter number, k = data point number

- redundancies
- virtual number of statistically independent data points with equal information content

$$R_p = \sum_s B_{sp}^i$$

s = spectrum number, p = parameter number




Fig. 1: File cut example of measurement with number density error

Analysis procedure and results

As a first application previously done measurements are re-analyzed with special focus on non-Voigt effects. The following analysis is based on a speed-dependent Voigt lineshape and does not account for temperature dependencies. A complete analysis is to be done in the near future.

1. Generation of a linestrength reference

- multispectrum fit of 9 room temperature pure water spectra
- Voigt-fit of position, linestrength, width parameter

2. Fit of column density and temperature for water/air-mixture room-temperature spectra

- single spectrum fit of mixture-spectra
- Voigt fit of position, effective linestrength, width
- fit of number density and temperature -> effective linestrengths match linestrength reference on average

3. Fit of all room-temperature spectra

- multispectrum fit of 9 pure and 16 mixture spectra (corrected column density and temperature)
- speed-dependent Voigt (SDV) lineshape for lines with a wide range of opacities - including opaque lines

Results

- opaque and non-opaque lines fitted simultaneously
- case 1: Voigt, case 2: speed-dependent Voigt (SDV)
- Voigt shows W-shaped residuals for non-opaque lines and line wing residuals for opaque lines
- speed-dependent Voigt improves residuals significantly




Fig. 3: Comparison between Voigt and SDV fitting residuals (residua are scaled for easier comparison)

- fitted broadening parameters systematically greater when fitting with SDV in comparison to a Voigt fit of lines with opacities < 4 [1] -> opaque lines are modeled too narrow with Voigt [3]
- influence of narrowing is greater when the broadening parameter is lower, i.e. when the J-quantum number is higher




Fig. 4: Relative differences of broadening parameters from SDV-analysis and HITRAN 2012




Fig. 5: Comparison of narrowing and broadening parameters

Conclusion

- a new multispectrum tool has been developed and tested against two independent tools
- innovations like an automatic mode and quality control mechanisms have been implemented
- widely used line profile models are implemented
- sophisticated ILS-modelling included
- as a first test water absorption spectra were re-analyzed
- Voigt lineshape not sufficient
- opaque and non-opaque lines have to be fitted simultaneously to prevent too narrow broadening parameters
- the narrowing parameter shows a clear trend with the broadening parameter
- complete analysis (T-dependencies, shifts, etc.) to be done

References

- [1] M. Birk, G. Wagner, Journal of Quantitative Spectroscopy and Radiative Transfer **113**, 889 (2012)
- [2] L. H. Coudert et. al., Journal of Molecular Spectroscopy **251**, 339 (2008)
- [3] G. Wagner, M. Birk, S. A. Clough, „Undiscovered errors of Voigt profile beyond tiny w-shaped residuals“, International Symposium on Molecular Spectroscopy – 68th Meeting, Columbus (OH), June 17th, 2013
- [4] G. Wagner, M. Birk, private communication, 06/2013
- [5] G. Wagner, M. Birk, F. Schreier, Journal of Geophysical Research **107**, ACH 10-1 - 10-18 (2002)
- [6] C. D. Boone, K. A. Walker, P. F. Bernath, Journal of Quantitative Spectroscopy and Radiative Transfer **105**, 525 (2007)
- [7] R. Ciurylo, J. Szudy, Journal of Quantitative Spectroscopy and Radiative Transfer **57**, 411 (1997)
- [8] M. Schneider et. al., Journal of Quantitative Spectroscopy and Radiative Transfer **112**, 465 (2011)
- [9] F. Hase, T. Blumenstock, C. Paton-Walsh, Applied Optics **38**, 3417 (1999)
- [10] J. Loos, M. Birk, G. Wagner, „Novel Design of MCT-Detector Optics with Enhanced Performance“, HRMS 2011, Dijon

Experiment

Experimental conditions

- Mar '03 - Jun '05
- Bruker IFS 120 HR
- improved MCT-detector [10]
- $p_{\text{total}} = 50 - 1000$ mb
- $p_{\text{H}_2\text{O}} = 0.001 - 5$ mb
- $l = 16 - 8500$ cm
- $T = 240 - 316$ K

Setup

- 26 cm double-jacket short cell
- double-jacket White-type multireflection cell with absorption path up to 85 m
- 800 l vessel for H₂O/air-mixture preparation
- Pt100 temperature sensors, mks baratron pressure gauges, mks constant pressure valve
- short cell: sealed off pure H₂O @ room temperature
- multireflection cell: pure H₂O flow @ room temperature
- multireflection cell: pressure broadened H₂O flow @ 240 - 316 K




Fig. 2: Experimental setup for water/air-mixture flow measurements

A more detailed description can be found in [1].

Acknowledgement: This work has been performed within the framework of DFG project „Improving spectroscopic data of H₂O and CH₄ for application on remote atmospheric measurements in the infrared“ under contract number BI 834/5-1.