

DLR-AG Facilities & Research Plans

Timo Frett¹, Jörn Rittweger¹, Edwin Mulder¹, Guido Petrat¹, Rupert Gerzer¹

(1) Institute of Aerospace Medicine, German Aerospace Center DLR, Linder Hoehe, 51147 Cologne, Germany

International Workshop on Research and Operational Considerations for Artificial Gravity Countermeasures

NASA Ames Research Center, Moffett Field, CA


20.02.2014

A large, curved image of the Earth from space occupies the bottom right portion of the slide. It shows a view of the Earth's surface with blue oceans, green landmasses, and white clouds. The curvature of the planet is clearly visible, and the text "Knowledge for Tomorrow" is overlaid on the right side of this image.

Knowledge for Tomorrow

AG Facilities at DLR – Overview and status

- Long Arm Human Centrifuge:



Closed in 2013 after more than 60 years of service!




AG Facilities at DLR – Overview and status

- ESA Short Arm Human Centrifuge:


Built	2008
Max. radius at outer perimeter	280 cm
Max. centrifugal acceleration	4.5 g
Max. g-alteration rate	0.325 g/s
Max. revolution of centrifuge rotor	39 rpm
Number and type of nacelles	2 beds 2 seats
Drive system	Three-phase motor
Power	75 kW
Max. payload	550 kg






Research facility :envifuge – *technical data*

Radius	380 cm
Max. acceleration	6g @ 370 cm
Max. onset rate (from 0g to 6g)	0.133 rad/s ² (≈ 0,2g/s)
Drive system	Torque drive, water-cooled
Power	30 kW
Max. angular velocity	4 rad/s (≈ 38 rpm)
Sound emission	< 40dB
Subjects height	150 – 210 cm
Max. subject weight	150 kg
Max. payload	800 kg
Nacelle translation speed	10 mm/s
No balancing needed. Extreme smooth running.	




© DLR



Research facility :envifuge – configurations



FLEXIBLE NACELLE

Adjustable nacelle for various configuration and positions.






Research facility :envifuge – nacelles

- Traveling distance of nacelles:
 - Main arm: 2,5 m
 - Secondary arm: 0,25-0,90 m (depending on subject height)
- Extremely short radius set-up possible:
 - Heart @ center of rotation
 - Head @ center of rotation
- Highly configurable:
 - Various mounting points
 - Auxiliary arms completely exchangeable
 - Auxiliary arms can be tilted



Research facility :envifuge – *robot arm*

- Axis:
 - 6 degrees of Freedom
- Speed
 - x, y: 30mm/s
 - z: 10mm/s
 - α, β, χ : 5°/s
- Payload
 - max. 1kg
- Position accuracy
 - +/- 1mm
- Contact force
 - Adjustable. Depends on probe geometry and safety settings.
- Interface
 - Active Contact Flange
 - Powerball



First Test Run – Jumps and Squats (1 g @ feet-level)



Future applications for space research on :envifuge

- Complex physical exercises under AG
 - Loading of bones and muscles
 - Neuro-vestibular effects (motion sickness)
 - Cardiovascular response with respect to the acceleration-gradient
- Comparison between Long arm and short arm centrifuge
- Short radius set up for realistic simulation of in-flight SRC's
 - Radius below 1.5m
 - Head position close to center of centrifugation
 - Multiple exercises (whole body training)
- Ultrasound examination of inner organs (e.g. heart) during centrifugation
- Both SRC's are available for National centrifuge program (NZP) and ESA Ground based facility program (CORA)
- International collaboration is very welcome

From a research- to a training-tool




For further information:
German Aerospace Center (DLR)
Institute of Aerospace Medicine
Linder Hoehe
51147 Cologne

Timo Frett
+49 2203 601 3458
timo.frett@dlr.de


Guido Petrat
+49 2203 601 3355
guido.petrat@dlr.de




Standard Nacelle




- The standard nacelle is designed for standard centrifuge runs. It can handle subjects with a body height from 1,50m to 2,10m and a body weight of 150kg at lateral or supine position.
- Additional windshields can be mounted to avoid unwanted air-stream.
- Instead of the force plate, a vibration plate can be mounted.



Jump platform




- Alternatively to the standard nacelle there is a lightweight sled available for performing squats and jumps.
- At the foot end there are two triaxial force plates which measures ground reaction force individually between left and right foot.




Motion Capturing System

- It consists of six cameras which covers almost 50% of the body. This system is mounted on a sled which can be moved along the lattice arm to position the MCS to capture the movement at the respective region of interest.


- VICON System
- 6 Bonita Cameras
- Body Builder Software
- Gigaswitch




Robotic arm



- The robotic arm is attached to a sled which can be positioned along the lattice arm. It consists of a mounting with a Powerball at the end.
- The Powerball has a clamping mechanism which can be equipped e.g. an ultrasound probe. This robotic arm is controlled from a special console in the control room.



Ergometer nacelle



- Lode Ergometer
- Remote Control
- Feedback

